

Goat Anti-Mouse IgG antibody, pre-adsorbed (PE)

Cat. No. GTX04199-08

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, ELISA
Reactivity	Mouse

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	≤ 1 µg/mL

Not tested in other applications.

Product Note

Reacts with the heavy and light chains of mouse IgG1, IgG2a, IgG2b, IgG2c, and IgG3 and with the light chains of mouse IgM and IgA. Pre-adsorbed with Human, rat, hamster, goat, sheep, rabbit, chicken, guinea pig, horse, and bovine serum proteins (SP). May react with immunoglobulins from other species

Properties

Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with mouse IgG
Purification	Purified by antigen-affinity chromatography
Conjugation	Phycoerythrin (PE) Wavelength


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04199-08 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04199-08 Goat Anti-Mouse IgG antibody, pre-adsorbed (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).