

Goat Anti-Mouse IgG1 (Heavy chain) antibody, F(ab')₂ fragment, pre-adsorbed**Cat. No. GTX04206**

Host	Goat
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	WB, ICC/IF, IHC-P, IHC-Fr, FCM, ELISA
Reactivity	Mouse

References (1)

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note

Pre-adsorbed with Mouse IgG2a, IgG2b, IgG3, IgM, and IgA; human immunoglobulins and pooled sera. May react with immunoglobulins from other species.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pepsin digest of Goat Anti-Mouse IgG1.
Purification	Purified IgG
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 16 Page 1 of 2


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04206 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04206 Goat Anti-Mouse IgG1 (Heavy chain) antibody, F(ab')₂ fragment, pre-adsorbed followed by Swine Anti-Goat IgG(H+L) antibody, pre-adsorbed (HRP).


GTX04206 FCM Image

FACS analysis of human peripheral blood lymphocytes using Mouse Anti-Human CD3 antibody followed by GTX04206 Goat Anti-Mouse IgG1 (Heavy chain) antibody, F(ab')₂ fragment, pre-adsorbed and Swine Anti-Goat IgG antibody, pre-adsorbed (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).