

Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (AP)

Cat. No. GTX04208-03

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Mouse

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:2000-1:4000

Not tested in other applications.

Product Note

Reacts with the heavy chain of C57BL/6 mouse IgG2c. Pre-adsorbed with Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA; human, bovine, and horse serum proteins (SP). May react with immunoglobulins from other species

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with mouse IgG2c
Purification	Purified by antigen-affinity chromatography
Conjugation	Alkaline Phosphatase (AP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04208-03 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04208-03 Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).