

## Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (Biotin)

Cat. No. GTX04209-02

|              |            |
|--------------|------------|
| Host         | Goat       |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | ELISA      |
| Reactivity   | Mouse      |

Package  
1 mg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ELISA              | 1:5000-1:20000       |

Not tested in other applications.

## Product Note

Reacts with the heavy chain of C57BL/6 mouse IgG2c. Pre-adsorbed with Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. May react with immunoglobulins from other species

## Properties

|               |                                                                                         |
|---------------|-----------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                  |
| Buffer        | PBS                                                                                     |
| Preservative  | 0.09% sodium azide                                                                      |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)              |
| Immunogen     | Pooled antisera from goats hyperimmunized with mouse IgG2c                              |
| Purification  | Purified by antigen-affinity chromatography                                             |
| Conjugation   | Biotin                                                                                  |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx04209-02 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04209-02 Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).