

Goat Anti-Mouse IgE (Heavy chain) antibody, pre-adsorbed (FITC)

Cat. No. GTX04222-06

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, FCM, IP, ELISA, IHC
Reactivity	Mouse

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:100-1:400
IHC	Assay dependent

Not tested in other applications.

Product Note Pre-adsorbed with Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. May react with immunoglobulins from other species.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pooled antisera from goats hyperimmunized with mouse IgE
Purification	Purified by antigen-affinity chromatography
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength


For full product information, images and publications, please visit our [website](#).

Date 2026 / 02 / 01 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04222-06 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04222-06 Goat Anti-Mouse IgE (Heavy chain) antibody, pre-adsorbed (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).