

## Rat Anti-Mouse IgD antibody [11-26] (Low endotoxin, azide free)

Cat. No. GTX04223-00

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Application	IHC-P, IHC-Fr, FACS, IP, Activation, PLA
Reactivity	Mouse

Package  
500 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FACS	Assay dependent
IP	Assay dependent
Activation	Assay dependent
PLA	Assay dependent

Not tested in other applications.

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG2a
Endotoxin	< 0.05 EU/µg (determined by the LAL assay)
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).