

Rat Anti-Mouse IgG2a antibody [SB84a] (AP)

Cat. No. GTX04227-03

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-Fr, FCM, ELISA, Multiplexing
Reactivity	Mouse

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	1:1000-1:4000
Multiplexing	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IgG2a hybridoma
Purification	Purified IgG1
Conjugation	Alkaline Phosphatase (AP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04227-03 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04227-03 Rat Anti-Mouse IgG2a antibody [SB84a] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).