

# Rat Anti-Mouse IgA antibody [11-44-2] (AP)

**Cat. No. GTX04228-03**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	WB, IHC-P, IHC-Fr, FCM, ELISA
<b>Reactivity</b>	Mouse

**Package**  
1 ml

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	1:1000-1:2000

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	50mM Tris, 1mM MgCl <sub>2</sub> , 50% Glycerol
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Purification</b>	Purified IgG1
<b>Conjugation</b>	Alkaline Phosphatase (AP)

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx04228-03 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04228-03 Rat Anti-Mouse IgA antibody [11-44-2] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).