

Rat Anti-Mouse IgA antibody [11-44-2] (APC)

Cat. No. GTX04228-07

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IHC-Fr, FCM, ELISA
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	$\leq 0.2 \mu\text{g}/10^6\text{cells}$
ELISA	$\leq 1 \mu\text{g/mL}$

Note : The suggested use of these reagents is in a final volume of 100µl.

Not tested in other applications.

Properties


Form	Liquid
Buffer	PBS, a stabilizer
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Allophycocyanin (APC) Wavelength

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04228-07 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04228-07 Rat Anti-Mouse IgA antibody [11-44-2] (APC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).