

## Rat Anti-Mouse lambda light chain antibody [JC5-1] (HRP)

Cat. No. GTX04229-01

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	WB, IHC-P, IHC-Fr, FCM, IP, ELISA, Activation, PLA
<b>Reactivity</b>	Mouse

Package

1 ml

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:4000-1:8000
Activation	Assay dependent
PLA	Assay dependent

**Note : Use it on non-reducing condition.**

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 50% glycerol
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Purification</b>	Purified IgG2b
<b>Conjugation</b>	Horseradish peroxidase(HRP)




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**  
 Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



**GTX04229-01 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04229-01 Rat Anti-Mouse lambda light chain antibody [JC5-1] (HRP).



For full product information, images and publications, please visit our [website](#).