

## Rat Anti-Mouse IgG2b antibody [SB74g] (HRP)

**Cat. No. GTX04231-01**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	WB, ICC/IF, FCM, ELISA
<b>Reactivity</b>	Mouse

**Package**  
1 ml

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent
ELISA	1:4000-1:8000

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 50% glycerol
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Mouse IgG2b hybridoma
<b>Purification</b>	Purified IgG2b
<b>Conjugation</b>	Horseradish peroxidase(HRP)


**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX04231-01 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04231-01 Rat Anti-Mouse IgG2b antibody [SB74g] (HRP).



For full product information, images and publications, please visit our [website](#).