

Rat Anti-Mouse IgG3 antibody [SB76b] (AP)

Cat. No. GTX04232-03

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM, ELISA
Reactivity	Mouse

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
ELISA	1:1000-1:2000

Not tested in other applications.

Properties

Form	Liquid
Buffer	50mM Tris, 1mM MgCl ₂ , 50% Glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IgG3 hybridoma
Purification	Purified IgG2a
Conjugation	Alkaline Phosphatase (AP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04232-03 ELISA Image

ELISA analysis of purified mouse immunoglobulins using serially diluted GTx04232-03 Rat Anti-Mouse IgG3 antibody [SB76b] (AP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).