

Goat Anti-Human IgG (Fab) antibody, pre-adsorbed (HRP)

Cat. No. GTX04236-01

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:2000-1:4000

Not tested in other applications.

Product Note

Reacts with the Fab region of human IgG, IgM, and IgA. Cross adsorption with human IgG Fc, may react with immunoglobulins from other species.

Properties

Form	Liquid
Buffer	PBS, 50% glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG, IgM, and IgA
Purification	Purified by antigen-affinity chromatography From serum
Conjugation	Horseradish peroxidase(HRP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04236-01 ELISA Image

ELISA analysis of purified human IgG Fab and IgG Fc proteins using serially diluted GTX04236-01 Goat Anti-Human IgG (Fab) antibody, pre-adsorbed (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).