

Influenza A virus H1N1 HA (Hemagglutinin) antibody [InA134]

Cat. No. GTX04240

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Influenza A virus (H1N1)

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX04239, Detection : GTX04240

Not tested in other applications.

Product Note Pair detects H1 of Influenza A as well as recombinant H1 with high specificity and does not interact with haemagglutinin H3.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09 % sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Influenza A/NewCaledonia/20/99 (strain H1N1) derived from allantoic fluid of 10 days old embryonated eggs
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES


GTX04240 ELISA Image

Sandwich ELISA analysis of viral lysates derived from different strains of Influenza A virus (i.e., Influenza A/New Caledonia/20/99(H1N1), Influenza A/Shangdong/9/93(H3N2), and mixture of Influenza B strains (B/Qingdao/102/91, B/Tokio/53/99, and B/Victoria/504/00)) using GTX04240 Influenza A virus H1N1 HA (Hemagglutinin) antibody [InA134].

Capture : GTX04239(1 µg/well)

Detection : GTX04240 (0.2 µg/well)



For full product information, images and publications, please visit our [website](https://www.genetex.com).