

AIRE antibody

Cat. No. GTX04257

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 58 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% sucrose
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human AIRE.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04257 WB Image


WB analysis of human lymphoblastoid cell lysates using GTX04257 AIRE antibody.
Dilution : 1:1000


GTX04257 IHC-P Image

IHC-P analysis of human kidney tissue using GTX04257 AIRE antibody.
Dilution : 4-8 µg/ml


GTX04257 WB Image

WB analysis of mouse skeletal muscle tissue lysates using GTX04257 AIRE antibody.
Dilution : 1 µg/mL


GTX04257 WB Image

WB analysis of human liver tissue lysates using GTX04257 AIRE antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).