

Hemoglobin beta S peptide

Cat. No. GTX04268-PEP

Species Human

Package

50 µg

Applications

Application Note

Hemoglobin beta S peptide can be used as a control with the GTX04268 Hemoglobin beta S antibody [23E5.H6.G6.C1.H7.F7.G9.F6]. Control peptide should be used at 1.0 µg per 1.0 µl of antiserum in per assay.

Properties

Form Lyophilized powder

Buffer Batch dependent

Storage Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Purity > 95%

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).