

Goat Anti-Rabbit IgG F(ab')₂ antibody (HRP)

Cat. No. GTX04274-01

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, IHC
Reactivity	Rabbit

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ELISA	1:10000-1:50000
IHC	1:500-1:2500

Not tested in other applications.

Product Note No reaction was observed against Rabbit IgG Fc.

Properties

Form	Liquid
Buffer	10mM Sodium Phosphate, 150mM NaCl, 1% BSA
Preservative	0.01% Gentamicin Sulfate
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit IgG F(ab') ₂ fragment
Purification	IgG fraction This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorptions to remove any unwanted reactivities.
Conjugation	Horseradish peroxidase(HRP)




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX04274-01 WB Image

WB analysis of various sample lysates using GTX04274-01 Goat Anti-Rabbit IgG F(ab')₂ antibody (HRP).

Lane 1 : Rabbit IgG whole molecule

Lane 2 : Rabbit IgG F(ab) Fragment

Lane 3 : Rabbit IgG Fc Fragment

Lane 4 : Rabbit IgM Whole Molecule

Lane 5 : Normal Rabbit Serum

Loading : 50 ng

Dilution : 1:10,000



For full product information, images and publications, please visit our [website](https://www.genetex.com).