

Glycated Hemoglobin A1c (HbA1c) antibody [75C9]

Cat. No. GTX04309

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Not tested in other applications.	

Product Note This antibody does not cross-react with non-glycated hemoglobin.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09 % sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Glycated (in vitro) synthetic peptide corresponding to the N-terminal end of beta-globin.
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04309 ELISA Image

ELISA analysis of native HbA1c or HbA1₀ proteins (100 ng/well) using GTX04309 Glycated Hemoglobin A1c (HbA1c) antibody [75C9].


GTX04309 ELISA Image

ELISA analysis of native HbA1c using GTX04309 Glycated Hemoglobin A1c (HbA1c) antibody [75C9] labeled with stable Eu³⁺ chelate (0.4 µg/well) as a detection antibody.

Capture antibody : GTX04311 Hemoglobin alpha antibody [Hb6] (1 µg/well).



For full product information, images and publications, please visit our [website](https://www.genetex.com).