

SMOC1 antibody

Cat. No. GTX04320

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:5000
ICC/IF	1:50-1:200
IHC-P	1:20-1:200

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.03% Proclin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Human SPARC-related modular calcium-binding protein 1 protein (277-382AA).
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04320 WB Image


WB analysis of mouse skeletal muscle tissue lysate using GTX04320 SMOC1 antibody.



GTX04320 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX04320 SMOC1 antibody.


Dilution : 1:100



GTX04320 IHC-P Image

IHC-P analysis of human colon cancer tissue using GTX04320 SMOC1 antibody.

Dilution : 1:100



GTX04320 IHC-P Image

IHC-P analysis of human gastric cancer tissue using GTX04320 SMOC1 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).