

Mouse Cardiac Troponin I protein, His tag

Cat. No. GTX04344-pro
Species Mouse

Package
50 µg

Properties
Form Lyophilized powder

Buffer Reconstitute with 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex. Lyophilized from 20mM Tris pH8.0, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL

Preservative 5% Trehalose and Proclin300

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Region/Sequence Full-length with N-terminal His Tag; Met1~Gly211 (NP_033432.1)

Expression System E. coli

Purity > 97%


Endotoxin <1.0 EU/1µg (determined by the LAL method)

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04344-pro Image

SDS-PAGE analysis of GTX04344-pro Mouse Cardiac Troponin I protein, His tag.



For full product information, images and publications, please visit our [website](https://www.genetex.com).