

Androgen Receptor antibody [MSVA-367R] HistoMAX™

Cat. No. GTX04368

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (1)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:100-1:200
-------	-------------

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to residues (within aa1-100) of Androgen Receptor.
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 09 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx04368 IHC-P Image

IHC-P analysis of human prostate tissue using GTx04368 Androgen Receptor antibody [MSVA-367R] HistoMAX™.


A strong nuclear androgen receptor immunostaining is seen in stromal and epithelial cells of the prostate.



GTx04368 IHC-P Image

IHC-P analysis of human colonic mucosa tissue using GTx04368 Androgen Receptor antibody [MSVA-367R] HistoMAX™.

Complete absence of androgen receptor staining in normal colon mucosa.



GTx04368 IHC-P Image

IHC-P analysis of human prostatic adenocarcinoma (PAC) tissue with (Gleason 4+4=8) tissue using GTx04368 Androgen Receptor antibody [MSVA-367R] HistoMAX™.

A strong nuclear AR immunostaining is seen in all cells of a prostatic adenocarcinoma (Gleason 4+4=8).



For full product information, images and publications, please visit our [website](#).