

TMEM16A antibody [MSVA-201M] HistoMAX™

Cat. No. GTX04369

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	IHC-P
Reactivity	Human

References (2)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution Recommended dilution

IHC-P 1:100-1:200

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer. Based on the testing result, the pH 9.0 showed a strongest staining pattern. The pH 7.8 is acceptable but lower pH results in a significant reduction of sensitivity.

This antibody has been validated by Autostainer including Agilent/Dako - Autostainer Link 48, Bond RX-Lecia, and Roche-Ventana Discovery ULTRA. The detailed steps by Autostainer please refer to the [protocol](#).

Not tested in other applications.

Product Note **Highly recommended for IHC-P in human tissues.**

Properties


Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human DOG-1 protein fragment (aa 2-101)



For full product information, images and publications, please visit our [website](#).


Date 2026 / 09 / 08 Page 1 of 2

Purification	Protein A/G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04369 IHC-P Image


IHC-P analysis of human gastrointestinal stromal tumor (GIST) tissue using GTX04369 TMEM16A antibody [MSVA-201M] HistoMAX™.

DOG1 Tumor periphery of a GIST showing intensive TMEM16A staining in all tumor cells.


GTX04369 IHC-P Image

IHC-P analysis of human seminal vesicle tissue using GTX04369 TMEM16A antibody [MSVA-201M] HistoMAX™.

Strong TMEM16A staining of apical membranes of glandular cells.


GTX04369 IHC-P Image

IHC-P analysis of human parotid gland tissue using GTX04369 TMEM16A antibody [MSVA-201M] HistoMAX™.

Strong TMEM16A staining of the apical membranes of secreting cells.



For full product information, images and publications, please visit our [website](https://www.genetex.com).