

Melan A antibody [MSVA-900M] HistoMAX™

Cat. No. GTX04383

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)

[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution Recommended dilution

IHC-P 1:100-1:200

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

This antibody has been validated by Autostainer including Agilent/Dako - Autostainer Link 48, Bond RX-Lecia, and Roche-Ventana Discovery ULTRA. The detailed steps by Autostainer please refer to the [protocol](#).

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues. This antibody is fully specific for Melan A. For Melan A+, we would recommend GTX04384 as a marker, which is for benign and malignant melanocytic tumors as well as steroid producing tissues.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human MART-1 protein




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 10 Page 1 of 2

Purification	Protein A/G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx04383 IHC-P Image

IHC-P analysis of human malignant melanoma from skin tissue using GTx04383 Melan A antibody [MSVA-900M] HistoMAX™.


Malignant melanoma with strong Melan A positivity in all tumor cells.



GTx04383 IHC-P Image

IHC-P analysis of human kidney tissue using GTx04383 Melan A antibody [MSVA-900M] HistoMAX™.

Melan A staining is absent in the kidney.



GTx04383 IHC-P Image

IHC-P analysis of human skin tissue using GTx04383 Melan A antibody [MSVA-900M] HistoMAX™.

A strong Melan A immunostaining is seen in melanocytes.



For full product information, images and publications, please visit our [website](#).