

Cytokeratin 13 antibody [MSVA-613M] HistoMAX™

Cat. No. GTX04385

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

References (1)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:100-1:200
-------	-------------

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer. Based on testing result, pH 9.0 showed a strongest staining pattern. However, pH 7.8 results in only a slight reduction of the staining intensity as compared to pH 9.0. Therefore, the pH 7.8 was an optimal for manual staining because of the better tissue preservation at pH 7.8 than at pH 9.0.

This antibody has been validated by Autostainer including Agilent/Dako - Autostainer Link 48, Bond RX-Lecia, and Roche-Ventana Discovery ULTRA. The detailed steps by Autostainer please refer to the [protocol](#).

Not tested in other applications.

Product Note	Highly recommended for IHC-P in human tissues.
--------------	--

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	esophageal keratins of rabbit origin



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 13 Page 1 of 2


Purification	Protein A/G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx04385 IHC-P Image

IHC-P analysis of human lip tissue using GTx04385 Cytokeratin 13 antibody [MSVA-613M] HistoMAX™. Strong cyto­keratin 13 staining in all suprabasal cell layers of the non keratinizing squamous epithelium of the lip. The basal cell layer is strictly cyto­keratin 13 negative.



GTx04385 IHC-P Image

IHC-P analysis of human skin tissue using GTx04385 Cytokeratin 13 antibody [MSVA-613M] HistoMAX™. Absence of cyto­keratin 13 staining in the keratinizing squamous epithelium of the skin.



GTx04385 IHC-P Image

IHC-P analysis of human oral cavity squamous cell carcinoma (OSCC) tissue using GTx04385 Cytokeratin 13 antibody [MSVA-613M] HistoMAX™. Strong cyto­keratin 13 immunostaining in a squamous cell carcinoma of the oral cavity.



For full product information, images and publications, please visit our [website](https://www.genetex.com).