

Cytokeratin 15 antibody [MSVA-615M] HistoMAX™

Cat. No. GTX04411

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-1:100

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer. Based on the testing result, the pH 9.0 showed a strongest staining pattern. The pH 7.8 is acceptable but lower pH results in a significant reduction of sensitivity.

Not tested in other applications.

Product Note **Highly recommended for IHC-P in human tissues.**

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human KRT15 protein
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 13 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04411 IHC-P Image


IHC-P analysis of human esophagus tissue using GTx04411 Cytokeratin 15 antibody [MSVA-615M] HistoMAX™.

A strong cytokeratin 15 immunostaining is always found in all cell layers of all non keratinizing and keratinizing squamous epithelium.


GTx04411 IHC-P Image

IHC-P analysis of human tonsil tissue using GTx04411 Cytokeratin 15 antibody [MSVA-615M] HistoMAX™.

A strong cytokeratin 15 immunostaining is found in all cell layers of all non keratinizing and keratinizing squamous epithelium except keratin layer.


GTx04411 IHC-P Image

IHC-P analysis of human tonsil tissue using GTx04411 Cytokeratin 15 antibody [MSVA-615M] HistoMAX™.

In tonsil crypts a strong cytokeratin 15 immunostaining is found in all cell layers of squamous epithelium.



For full product information, images and publications, please visit our [website](#).