

HLA-DRB1 antibody [MSVA-478R] HistoMAX™

Cat. No. GTX04457

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (2)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:100-1:200
-------	-------------

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human HLA-DRB1 fragment
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 11 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx04457 IHC-P Image

IHC-P analysis of human lung tissue using GTx04457 HLA-DRB1 antibody [MSVA-478R] HistoMAX™. Strong HLA-DRB1 staining in alveolar macrophages HLA-DRB1 immunohistochemistry.



GTx04457 IHC-P Image

IHC-P analysis of human testicular seminoma tissue using GTx04457 HLA-DRB1 antibody [MSVA-478R] HistoMAX™. HLA-DRB1 negative seminoma showing HLA-DRB1 staining of inflammatory cells and of small vessels HLA-DRB1 immunohistochemistry.



GTx04457 IHC-P Image

IHC-P analysis of human liver tissue using GTx04457 HLA-DRB1 antibody [MSVA-478R] HistoMAX™. Intense HLA-DRB1 staining of inflammatory cells and of Kupffer cells HLA-DRB1 immunohistochemistry.



For full product information, images and publications, please visit our [website](https://www.genetex.com).