

UPK3B antibody [MSVA-736M] HistoMAX™

Cat. No. GTX04476

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (1)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:100-1:200
-------	-------------

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

This antibody has been validated by Autostainer including Agilent/Dako - Autostainer Link 48, Bond RX-Lecia, and Roche-Ventana Discovery ULTRA. The detailed steps by Autostainer please refer to the [protocol](#).

Not tested in other applications.

Product Note **Highly recommended for IHC-P in human tissues.**

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human UPK3B fragment
Purification	Protein A/G purified



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2


Conjugation

Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04476 IHC-P Image


IHC-P analysis of human epithelioid mesothelioma tissue using GTx04476 UPK3B antibody [MSVA-736M] HistoMAX™.

Epithelioid mesothelioma with membranous UPK3B staining of variable intensity in most tumor cells.


GTx04476 IHC-P Image

IHC-P analysis of human peritoneal mesothelium tissue using GTx04476 UPK3B antibody [MSVA-736M] HistoMAX™.

Distinct Upk3b immunostaining of apical membranes of peritoneal mesothelial cells.


GTx04476 IHC-P Image

IHC-P analysis of human urinary bladder tissue using GTx04476 UPK3B antibody [MSVA-736M] HistoMAX™.

Urothelium with distinct UPK3B immunostaining of the apical membranes of umbrella cells.



For full product information, images and publications, please visit our [website](https://www.genetex.com).