

Rabbit anti-human IgA antibody [MSVA-700R] HistoMAX™

Cat. No. GTX04490

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

IHC-P

1:100-1:200

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human IGHA1 and IGHA2 proteins
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 12 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx04490 IHC-P Image

IHC-P analysis of human mucosa from duodenum tissue using GTx04490 Rabbit anti-human IgA antibody [MSVA-700R] HistoMAX™.


Human IgA positive plasma cells are abundant in the duodenum mucosa.



GTx04490 IHC-P Image

IHC-P analysis of human tonsil tissue using GTx04490 Rabbit anti-human IgA antibody [MSVA-700R] HistoMAX™.

Human IgA positive plasma cells are numerous in the tonsil.



GTx04490 IHC-P Image

IHC-P analysis of human squamous cell carcinoma (SCC) from skin tissue using GTx04490 Rabbit anti-human IgA antibody [MSVA-700R] HistoMAX™.

Squamous cell carcinoma containing many human IgA positive plasma cells in its stroma.



For full product information, images and publications, please visit our [website](#).