

Rabbit anti-human IgD antibody [MSVA-701R] HistoMAX™

Cat. No. GTX04491

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)

[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Note : Manual staining : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IGHD fragment
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 12 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx04491 IHC-P Image

IHC-P analysis of human tonsil tissue using GTx04491 Rabbit anti-human IgD antibody [MSVA-701R] HistoMAX™.


Strong human IgD staining in a fraction of B lymphocytes especially in the mantle zone.



GTx04491 IHC-P Image

IHC-P analysis of human diffuse large B-cell lymphoma (DLBCL) from lymph node tissue using GTx04491 Rabbit anti-human IgD antibody [MSVA-701R] HistoMAX™.

Diffuse large B cell lymphoma with intense human IgD immunostaining of tumor cells.



GTx04491 IHC-P Image

IHC-P analysis of human thymus tissue using GTx04491 Rabbit anti-human IgD antibody [MSVA-701R] HistoMAX™.

Only few cells show human IgD staining especially in the medulla.



For full product information, images and publications, please visit our [website](https://www.genetex.com).