

Human p53 protein, His tag

Cat. No. GTX04503-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

Applications

Observed MW (kDa)	57 kDa.
-------------------	---------

Properties

Form	Lyophilized powder
------	--------------------

Buffer	Reconstitute with ddH ₂ O to 0.1-0.2 mg/mL. Do not vortex. Lyophilized from 20mM Tris (pH8.0), 150mM NaCl, 0.01% SKL, 5% Trehalose.
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
---------	--

Region/Sequence	N-terminal His tagged full length human p53 protein (Met1~Asp393, Uniprot #P04637-1)
-----------------	--

Expression System	E. coli
-------------------	---------

Purity	>90%
--------	------


Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX04503-pro Image

SDS-PAGE analysis of GTX04503-pro Human p53 protein, His tag.



For full product information, images and publications, please visit our [website](#).