

Human H-Ras (G12V mutant) protein, His tag

Cat. No. GTX04505-pro

Species	Human
----------------	-------

Package	50 µg
----------------	-------

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Reconstitute with ddH ₂ O to 0.2 mg/mL. Do not vortex. Lyophilized from 10 mM Hepes (pH7.4), 500 mM NaCl, 5% trehalose.
---------------	--

Preservative	No preservatives
---------------------	------------------

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Region/Sequence	N-terminal His tagged full length human HRAS G12V protein (1-189aa, Uniprot #P01112-1)
------------------------	--

Expression System	E. coli
--------------------------	---------

Purity	>95% by SDS-PAGE
---------------	------------------

Conjugation	Unconjugated
--------------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).