

## Lipid A LPS antibody [26-5]

Cat. No. GTX04515

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	Immunoassay
Reactivity	Bacteria

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunoassay	1:10

Not tested in other applications.

**Product Note** This antibody recognizes Lipid A on the outer membrane of most gram-negative bacteria.

## Properties

Form	Liquid
Buffer	Filter-sterilized PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein G purified
Conjugation	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).