

Gasdermin D antibody [6D11]

Cat. No. GTX04538

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.25-0.5 µg/mL
ICC/IF	5 µg/mL
IHC-P	2-5 µg/mL
FCM	1-3 µg/1x10 ⁶ cells
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived human GSDMD recombinant protein (Position: M1-H484).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04538 FCM Image


FACS analysis of Jurkat cells using GTx04538 Gasdermin D antibody [6D11].

Blue : Primary antibody

Green : Isotype control

Red : Unlabelled sample


Dilution : 1 µg/1x10⁶ cells


GTx04538 IHC-P Image

IHC-P analysis of human esophageal squamous carcinoma tissue using GTx04538 Gasdermin D antibody [6D11].

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)


Dilution : 2 µg/mL


GTx04538 WB Image

WB analysis of various samples using GTx04538 Gasdermin D antibody [6D11].

Loading : 30 µg

Dilution : 0.5 µg/mL


GTx04538 IHC-P Image

IHC-P analysis of human cervical cancer tissue using GTx04538 Gasdermin D antibody [6D11].

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).