

SLC20A1 antibody

Cat. No. GTX04583

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.25-0.5 µg/mL
IHC-P	2-5 µg/mL
ELISA	0.1-0.5 µg/mL

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived human SLC20A1 recombinant protein (Position: A79-K494).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04583 IHC-P Image

IHC-P analysis of human lung adenocarcinoma tissue using GTx04583 SLC20A1 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)


Dilution : 2 µg/mL


GTx04583 IHC-P Image

IHC-P analysis of mouse cerebellum tissue using GTx04583 SLC20A1 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)


Dilution : 2 µg/mL


GTx04583 IHC-P Image

IHC-P analysis of mouse colon tissue using GTx04583 SLC20A1 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL


GTx04583 IHC-P Image

IHC-P analysis of rat kidney tissue using GTx04583 SLC20A1 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).