

Human LDL protein

Cat. No. GTX04613-pro

Species Human

Package
10 mg

PRODUCT

Summary The physical appearance of Human LDL protein is yellow to orange liquid. And human LDL protein can be further diluted with PBS containing 15% sucrose.

Properties

Form Liquid

Buffer Batch dependent

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence Native protein

Expression System Human plasma

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).