

SNF5 antibody [ARC53139]

Cat. No. GTX04623

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:5000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 44 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	0.05% proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SMARCB1/SNF5 (NP_003064.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04623 IHC-P Image

IHC-P analysis of human undifferentiated carcinoma of esophagus tissue using GTX04623 SNF5 antibody [ARC53139].

Antigen retrieval : 10 mM Tris/EDTA buffer pH 9.0


Dilution : 1:100


GTX04623 WB Image

WB analysis of various samples using GTX04623 SNF5 antibody [ARC53139].

Loading : 25 µg


Dilution : 1:10000


GTX04623 WB Image

WB analysis of HeLa cells with/without KD treatment using GTX04623 SNF5 antibody [ARC53139].

Loading : 25 µg

Dilution : 1:2000


GTX04623 IHC-P Image

IHC-P analysis of human epithelioid sarcoma (ini-1 deletion) tissue using GTX04623 SNF5 antibody [ARC53139].

Antigen retrieval : 10 mM Tris/EDTA buffer pH 9.0

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).