

## JNK antibody [DFI-13]

Cat. No. GTX04659

|              |                                                   |
|--------------|---------------------------------------------------|
| Host         | Rabbit                                            |
| Clonality    | Monoclonal                                        |
| Isotype      | IgG                                               |
| Applications | WB, ICC/IF, IP                                    |
| Reactivity   | Human, Mouse, Rat, Zebrafish, Bovine, Dog, Monkey |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000-1:2000        |
| ICC/IF             | 1:50-1:200           |
| IP                 | 1:50                 |

Not tested in other applications.

Calculated MW 42-47 kDa. ([Note](#))

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 150mM NaCl, 0.04-0.05% BSA, 50% Glycerol                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | A synthesized peptide derived from human JNK1 + JNK2 + JNK3                                                                                                                                                                |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx04659 WB Image

WB analysis of HeLa whole cell lysate using GTx04659 JNK antibody [DFI-13].



### GTx04659 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx04659 JNK antibody [DFI-13].

Green : Primary antibody

Blue : DAPI



For full product information, images and publications, please visit our [website](https://www.genetex.com).