

HSP22 antibody

Cat. No. GTX04663

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:100
IHC-P	1:100
IP	Assay dependent

Not tested in other applications.

Calculated MW 22 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% glycerol
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human HSP22
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx04663 WB Image

WB analysis of rat skeletal muscle tissue lysate using GTx04663 HSP22 antibody.

Loading : 15 µg

Dilution : 1:1000




GTx04663 IHC-P Image

IHC-P analysis of mouse backskin tissue using GTx04663 HSP22 antibody.

Localization : Epidermis positive, dermal staining.

Fixation: Bouin's Fixative Solution.

Dilution : 1:100



GTx04663 ICC/IF Image

ICC/IF analysis of cold 100% methanol-fixed HaCaT cells using GTx04663 HSP22 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).