

Human CD66e protein, Fc tag

Cat. No. GTX04664-pro

Species Human

Package

50 µg

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O to a concentration less than 100 µg/ml. Lyophilized from a 0.2 µm filtered solution of PBS pH 7.4.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Region/Sequence C-terminal Fc tagged, mature human CD66e protein (P06731, aa Lys35-Ala685)

Expression System Human cells

Purity >95% by SDS-PAGE

Endotoxin < 1 EU/µg (Determined by LAL assay)

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).