

Dengue virus type 2 Envelope protein

Cat. No. GTX04675-pro

Species Dengue virus 2

Package
100 µg

Properties

Form Liquid

Buffer PBS

Preservative 0.1% ProClin300

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length Dengue type 2 envelope protein

Expression System Insect cells

Purity >95% by SDS-PAGE

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).