

Helicobacter pylori CagA antibody [HP-1811]

Cat. No. GTX04677

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, EIA
Reactivity	Helicobacter pylori

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
EIA	Assay dependent

Not tested in other applications.

Product Note This antibody recognizes epitope localized between Asp562 and Gln795 in the GenBank FJ428215 amino acid sequence.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09 % sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragments of CA-gene of H. pylori.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).