

Leptospira interrogans LipL41 antibody [0412]

Cat. No. GTX04681

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ICC/IF, ELISA, Sandwich ELISA
Reactivity	Leptospira interrogans, Leptospira kirschneri

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200
Sandwich ELISA	Assay dependent

Note : Capture : GTX04682 / GTX04683 / GTX04687, Detection : GTX04681

Not tested in other applications.

Product Note

It reacts with: 4 Serovars of Leptospira interrogans + 2 Serovars of L. kirschneri and does not react with Leptospira biflexa (as expected) and mammalian cells.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).