

Horseradish Peroxidase antibody

Cat. No. GTX04689

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Horseradish

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:4000-1:8000
ELISA	1:4000-1:8000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified horseradish peroxidase
Purification	Purified by affinity chromatography From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04689 WB Image

WB analysis of purified horseradish peroxidase protein using GTX04689 Horseradish Peroxidase antibody.

Lane 1 : 0.2 µg of purified HRP

Lane 2 : 0.1 µg of purified HRP

Dilution : 1:8000



For full product information, images and publications, please visit our [website](#).