

CD37 antibody [MB-1]

Cat. No. GTX04698

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-Fr, FCM, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
FCM	1-4 µg/mL
IP	Assay dependent

Note : Do not need a permeabilization step.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	15 mM sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human lymphomatous lymph node cells
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX04698 FCM Image

FACS analysis of human CD37 positive lymphocytes separated from human peripheral whole blood cells using GTX04698 CD37 antibody [MB-1].


Red-filled : CD37 positive

Black-dashed : CD37 negative

Dilution : 0.2 µg/mL


GTX04698 FCM Image

FACS analysis of human lymphocytes using GTX04698 CD37 antibody [MB-1] (0.2 µg/mL) and GTX18229 CD19 antibody [LT19] (PE) (20 µL reagent / 100 µL of peripheral whole blood).


GTX04698 FCM Image

FACS analysis of human peripheral whole blood cells using GTX04698 CD37 antibody [MB-1].

Dilution : 0.2 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).