

CD79a antibody [4A10D12]

Cat. No. GTX04701

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:200-1:1000
FCM	1:200-1:400
ELISA	1:10000

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human CD79A (AA: extra(33-143)) expressed in E. Coli.
Purification	Purified IgG2b
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04701 FCM Image

FACS analysis of Raji cells using GTx04701 CD79a antibody [4A10D12].

Green : Primary antibody

Red : Negative control


GTx04701 WB Image


WB analysis of various samples using GTx04701 CD79a antibody [4A10D12].

Lane 1 : Raji whole cell lysate

Lane 2 : Ramos whole cell lysate

Lane 3 : K562 whole cell lysate

Lane 4 : Molt-4 whole cell lysate


GTx04701 WB Image

WB analysis of HEK293 cells with/without CD79a (aa 33-143)-hlgGfc overexpression using GTx04701 CD79a antibody [4A10D12].



For full product information, images and publications, please visit our [website](https://www.genetex.com).