

Cytokeratin 5/14 antibody [XM26 + LL002]

Cat. No. GTX04713

Host	Mouse
Clonality	Multiclonal
Isotype	IgG1, IgG3
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:25-1:150

Note : Recommendations for antigen retrieval : Citrate Buffer pH 6.0

Not tested in other applications.

Properties

Form	Liquid
Buffer	batch dependent
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	XM26: Recombinant protein corresponding to 103 amino acid portion of C-terminal region of the human cytokeratin 5 molecule_x000D_ LL002: Synthetic peptide of the extreme C-terminal of human keratin-14
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04713 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX04713 Cytokeratin 5/14 antibody [XM26 + LL002].
Antigen retrieval : Citrate Buffer pH 6.0



For full product information, images and publications, please visit our [website](#).