

NDUFS2 antibody

Cat. No. GTX04733

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.25 µg/mL
ICC/IF	5 µg/mL
IHC-P	2-5 µg/mL
FCM	1-3 µg/10 ⁶ cells
ELISA	0.1-0.5 µg/mL

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived human NDUFS2 recombinant protein (Position: M1-R463).
Purity	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX04733 IHC-P Image

IHC-P analysis of mouse heart tissue using GTX04733 NDUF52 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL


GTX04733 IHC-P Image


IHC-P analysis of human liver cancer tissue using GTX04733 NDUF52 antibody.

Red : Primary antibody

Blue : DAPI

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)


Dilution : 5 µg/mL


GTX04733 IHC-P Image

IHC-P analysis of human breast cancer tissue using GTX04733 NDUF52 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL


GTX04733 IHC-P Image

IHC-P analysis of human liver cancer tissue using GTX04733 NDUF52 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 2 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).