

Progesterone Receptor antibody [SP2]

Cat. No. GTX04737

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:25-1:100

Note : Recommendations for antigen retrieval : EDTA Buffer pH 8.0 at room temperature for 30-60 mins.

Not tested in other applications.

Properties

Form	Liquid
Buffer	batch dependent
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encoding human progesterone receptor 412-526 aa.
Purification	Unpurified From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04737 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX04737 Progesterone Receptor antibody [SP2].
Antigen retrieval : EDTA Buffer pH 8.0



For full product information, images and publications, please visit our [website](#).