

Cytomegalovirus antibody [DDG9 + CCH2]

Cat. No. GTX04739

Host	Mouse
Clonality	Multiclonal
Isotype	IgG1, IgG2a
Applications	IHC-P
Reactivity	Cytomegalovirus

References (1)
Package
6 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:10-1:50

Note : Recommendations for antigen retrieval : Citrate Buffer pH 6.0

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Cytomegalovirus infected cell lysate.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04739 IHC-P Image

IHC-P analysis of Cytomegalovirus infected human tissue using GTX04739 Cytomegalovirus antibody [DDG9 + CCH2].

Antigen retrieval : Citrate Buffer pH 6.0



For full product information, images and publications, please visit our [website](#).