

FBXO39 antibody, C-term

Cat. No. GTX04746

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

References (1)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	a KLH conjugated synthetic peptide between 385-414 amino acids from the C-terminal region of human FBXO39.
Purification	Protein A purified followed by peptide affinity purification
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx04746 WB Image

WB analysis of mouse kidney tissue lysate using GTx04746 FBXO39 antibody, C-term.

Loading : 35 µg



GTx04746 IHC-P Image

IHC-P analysis of human brain tissue using GTx04746 FBXO39 antibody, C-term.



For full product information, images and publications, please visit our [website](https://www.genetex.com).